



# Rotaid AED Cabinets

## Accessoires - Totem



English  
version

# Product sheet

*The Rotaaid AED Cabinets totem allows for a stand alone in- and outdoor installation of your defibrillator. The solid construction and high quality appearance make sure your defibrillator stands out and can be securely placed in practically all environments and situations.*

The modular set-up ensures an easy installation and the totem can easily be demounted when used for training and trade show purposes.

On the front, back and top of the totem, information and branding panels can be mounted. These panels are custom made, connecting to the customer's request and location, making sure the totem aligns perfectly with its surroundings.

The cabinet is safely mounted on the totem and the solid steel construction allows for installation in the most demanding environments. The double powder coated finish ensures a long-lasting appearance.

## Appliances

- Training and fares
- Stand alone installations
- Locations with no wall mounting possibilities
- Outdoor installations
- Renting and eventing

## Design

- Custom made design
- Full colour UV protected printing
- Allows for fixed floor mounting
- In combination with all Rotaaid cabinets

- ✚ For indoor and outdoor usage
- ✚ Solid and safe construction
- ✚ Storage box for transformer unit
- ✚ Custom made designs and panels



# Specifications

## Dimensions

### **Total**

Width: 700 mm

Height: 1615 mm

Weight: 40 kg

### **Footplate**

Width: 700 mm round

Thickness: 5 mm high

Material: Steel, powder-coated

### **Uprights / Stairs**

Width: 560 mm

Height: 1715 mm

Thickness pole: 60 mm

Material: Steel, powder-coated

### **Panel**

Width: 394 mm

Height: 1075 mm

Material: Dibond, 3 mm thick

## Modules

### **1. Sunroof**

The perfect protection to maximize shade.

### **2. Backboard** *(Custom made)*

Provide additional guidance and instructions.

### **3. Panel** *(Custom made)*

Promote your totem with a custom made sign.

### **4. Additional color**

Available in every color.

## Article numbers

### **9000**

Rotaid Totem with foot and front sign.

### **9000 + 302**

Rotaid Totem with foot, front sign and backboard.

### **9001**

Rotaid Totem with foot, front sign and sunroof.



**ROTAID™**  
AED CABINETS



### Article number - 9000

Rotaid Totem with foot and front sign.



### Article number - 9000+302

Rotaid Totem with foot, front sign and backboard.



### **Article number - 9001**

Rotaid Totem with foot, front sign and sunroof.