



CardiLink

CONNECTED
AED SAFETY
SOLUTIONS





Sudden Cardiac Arrest is one of the world's leading causes of death. Every minute without using a defibrillator the chance of survival decreases by seven to ten per cent. With a deployable Automated External Defibrillator, anyone can help saving lives.

CardiLink® state of the art connectivity solutions ensures 24/7 AED readiness for CardiAid AEDs.



Always Ready

Checks & Documents AED status



Automated QMS Compliance

Certified & Online documentation



Emergency Notifications

Sends a notification when moved & opened



Uptime Guarantee

Via pro-active servicing & SLA



GPS Geofence

Theft detection movement tracking



Fleet Management

Managed solutions, API integrated




CardiLink

Drahtlose Konnektivitätslösungen für CardiAid AEDs

CardiAid®
Automated External Defibrillator

Providing connected
safety solutions
as a Service



In today's world, customers expect a safety solution that works! And if you are the customer and user of the solution, why should YOU have to worry if the solution works? The provider of the solution should take care of this!

- Daily readiness updates
- SMS alerts when device is in motion



Get connected with as many AEDs as required.



Get advice from the experts on safe zones and AED placement.



Get the reports in your own IT system via API integration.



Get comfortable,
We take care.

- YOU need to walk around to check equipment?

That is the world upside down!

Let US remotely monitor this for you. Saving you time, money and organizational headache

- YOU need to call service?

That is the world upside down!

Let US inform you that an issue was detected at night which was fixed in the morning

- YOUR employees & customers unsafe while YOU wait for service?

That is the world upside down!

Let US fix it before rush hour so everyone is safe! No need for delays and unsafe situations

- YOU having to document everything?

That is the world upside down!

Let US document online and ensure your compliance

- YOU need to find this documentation in yet another IT system?

That is the world upside down!

Let US integrate the documentation and any ticket open and close with your current IT solution

- YOU need to worry about theft?

That is the world upside down!

Let US send notifications to your security guard with GPS coordinates when the AED goes beyond a pre-set geo fence

- YOU wouldn't know the AED was moved and then opened?

That is the world upside down!

Let US send notifications to your first aid team so they know where the AED went and that it is opened

Data transmitting
through secured and
encrypted networks

CardiAid®
Automated External Defibrillator



Your distributor in Switzerland

MSTRC Medical sàrl
Chemin de Préveyres 57e
CH-1132 Lully sur Morges (VD)
+41 22 575 39 90
www.mstrc-medical.ch | info@mstrc-medical.ch

