

Escape-Chair® VOLT

Powered stair climbing evacuation & rescue chair...



What's an Escape-Chair® VOLT?

An Escape-Chair® VOLT is the **POWERED** solution to evacuate & rescue people to safety in a seated position up and/or down a flight of stairs. A simple procedure can make the Escape-Chair® VOLT fully functional within a few seconds.

How does the Escape-Chair® VOLT work?

The Escape-Chair® VOLT is operated by a self-explanatory control panel with display on the operating handle. After unfolding the Escape-Chair® VOLT it is immediately ready for use. The single operator must push the operating handle down in the direction of the stairs.

Advantages of the Escape-Chair® VOLT:

- Powered support makes going up the stairs much easier for operator.
- Ideal for quick evacuation (4 speed descent).
- Lithium cell enough for up/down 40 floors.
- Single operator chair.
- Ready for use within seconds.
- Easy to use and operate (wireless).
- Equipped with arm- and footrest.
- 4 height adjustable operating handle.
- Battery indicator on display control panel.
- LED light for dark environment evacuation.
- CE certified.

Our advice:

Make sure you have an Escape- Chair® VOLT available in your organization so that you're prepared for any evacuation and/or emergency, for up and/or down evacuation.

"Failure to prepare for an emergency has its consequences."





Specifications	Escape-Chair® VOLT
Article number:	10050
Dimensions folded up (height x width x depth):	109 x 63 x 29 cm
Maximum load bearing:	180 kg
Maximum stair gradient:	40°
Weight:	±37,5 kg
Using time when full charged:	2 hours (± 40 floors)
Type of tracks:	Serrated tracks for superior grip
Unlocking function tracks:	✓
Battery life:	2 – 3 years
Extra batterypack:	Optional
Number of speeds:	4
Ready for use within seconds:	✓
Foot rest:	✓
Arm rests:	✓
LED lightning for dark environment:	✓
4 height adjustable operating handle:	✓
Flame & heat resistant dust cover:	Optional
Carrying handles front & back:	✓
10G crash test bracket for Ambulances / vehicles	✓
CE certified:	✓
Wall mounting bracket included	✓
Delivery incl. wall bracket:	✓
Warranty on non-wearing parts:	1 year (Battery 6 months)